


## Mixing Instructions

### Step 1 | Unpacking

- **Open packaging** and transfer contents to the sterile field using sterile technique.
- **Loosen** stopper on each 1 mL syringe by pulling plunger rod back slightly to break seal.

### Step 2 | Assembly

- **Remove the cap** from the prefilled mixing syringe containing the powder (**do not discard**).




- **Remove the cap** from the glass syringe.

- Using a female-female Luer lock connector, attach the glass syringe containing the setting solution to the mixing syringe containing the powder.

### Step 3 | Mixing

- Hold the syringes vertically with the glass syringe on top. **Inject entire contents of setting solution** into the mixing syringe containing the powder. **Pull back** powder syringe plunger rod to ensure contents are fully transferred.



- **Start timer** as soon as powder and setting solution make contact.
- Remove the solution syringe and Luer lock connector (**do not discard**). **Recap** the mixing syringe.
- Remove the plunger rod sleeve (**do not discard**) and use the integrated mixing handle to mix the solution and powder by pressing the handle in and out of the material from top to bottom, and rotating the handle if necessary.
- **Mix for 1 minute** to ensure complete mixing.



### Step 4 | Loading injection syringes

- **Reattach** the plunger rod sleeve over the integrated mixing handle. Remove tip cap and attach a Luer lock connector to the mixing syringe. **Expel residual air** from the mixing syringe until a small drop of Tactoset is ejected.
- Attach a 1 mL Luer lock syringe. **Hold and pull back** on the plunger rod of the 1 mL syringe to aid in transferring Tactoset material. Fill to 0.9 mL. **Do not overfill.**
- Fill and repeat as needed with remaining syringes.



## Timing Notes

*Start timer as soon as powder and setting solution make contact*

0:00

### Minute 0-1 | Mixing

- Mix powder and solution continuously for 1 minute

1:00

### Minutes 1-7 | Fill syringes & wait

- Fill the provided 1 mL syringes to 0.9 mL
- Wait to inject until at least minute 7

7:00

### After Minute 7 | Injection

- Inject Tactoset at desired timepoint after 7 minutes
- Continue checking material for paste-like consistency to ensure injectability
- Inject at surgeon's desired consistency
- Apply continuous pressure when injecting
- Inspect for extrusion and remove any excess material within 2 minutes following implantation

## Ordering Information

Item Code	Description
1000040	<b>Medium bone 11Ga delivery cannula set</b> including: 1 drillable outer cannula with removable drillable stylet; 1 inner cannula for end delivery; 1 inner cannula for side delivery; 1 push rod
4000050	<b>Small bone 15Ga delivery cannula set</b> including: 1 drillable outer cannula with removable drillable stylet; 1 push rod
6000041	<b>Tactoset material set</b> including: 1 prefilled mixing syringe containing Tacoset powder; 1 glass syringe containing Tactoset setting solution; 4 Luer lock syringes (1 mL); 1 female-female Luer lock connector

*Please note that both the Tactoset material kit item code and one of the delivery cannula item codes are needed for each procedure.*

For complete product information, please visit <http://www.anikaifu.com/>

Anika and Tactoset are registered trademarks of Anika Therapeutics, Inc.  
 ©2020 Anika Therapeutics, Inc. | 07/20 | AML 900-122/B

**Anika Therapeutics, Inc.**  
 32 Wiggins Ave, Bedford, MA  
 1-888-721-1600  
[www.anikatherapeutics.com](http://www.anikatherapeutics.com)