

The Synd-EZ kit provides a solution for damage to the syndesmosis joint. This simple yet effective device is designed to stabilize the ankle and can be used with or without a fracture plate.

## Versatile implant design

- Ideal for use with or without a fracture plate
- Washer included for use without a fracture plate\*
- Top hat is low profile and self centering on fracture plates



## Adjustable knotless construct

- Simplifies surgical procedure
- Offers confidence of flipped medial button from white/black toggle suture

## Counter traction suture & tension handle

- Integrated for predictably smooth tension of construct




\*Washer is only available in the Synd-EZ Titanium Kit 11223.

# Synd-EZ Kit *Syndesmosis Repair*

## SYND-EZ KIT COMPONENTS

1. Top hat
2. White/blue adjustment suture
3. Solid blue counter-traction suture
4. Button
5. Solid white passing suture
6. White/black toggle suture
7. 3.5mm cannulated drill bit
8. 2.4mm passing pin
9. 1.5mm K-wire
10. Black tension handle
11. Washer\*



## SYND-EZ ORDERING INFORMATION

Synd-EZ				
Part #	Width of Button (mm)	Length of Button (mm)	Description	
11223	3.2	12	Synd-EZ Titanium Kit (Ti)	
11224	3.2	12	Synd-EZ Stainless Steel Kit (SS)	

Synd-EZ Instrumentation (for use with a fracture plate)					
Part #	Description			Sterile/Non-sterile	Single-use/Reusable
11272	Synd-Ez Drill Guide			Non-sterile	Reusable

\* Washer is only available in the Synd-EZ Titanium Kit 11223.



Parcus Medical has joined **Anika**

**Parcus Medical, LLC.**

6423 Parkland Dr, Sarasota, FL

1-941-755-7965 · [www.parcusmedical.com](http://www.parcusmedical.com)

[ParcusCS@anika.com](mailto:ParcusCS@anika.com)

©2021 Parcus Medical, LLC. | QD 8042 Rev 3 | AML 900-227/A | 12/2021  
Anika and Parcus are registered trademarks of Anika Therapeutics and their affiliates

