

# FREEFIX<sup>®</sup>

ulnar shortening system



The freedom to perform a distal or proximal osteotomy  
based on clinical need

# A system that allows the distal interosseous oblique ligament to dictate the osteotomy position

## Shortening options to address specific pathology


### DISTAL OSTEOTOMY

metaphyseal location facilitates healing and maintains normal distal interosseous ligament tension

### PROXIMAL OSTEOTOMY

increases tension on distal interosseous ligament to enhance DRUJ stability

## Plate features




# Custom instrumentation for an efficient, reproducible procedure

Distal and proximal cutting guides allow blade access from either side of the plate

## DISTAL PLATE GUIDE

Transverse Cut




## PROXIMAL PLATE GUIDES

Oblique or Transverse Cuts



Streamlined procedure that maintains rotational stability





**skeletal dynamics**<sup>®</sup>  
UNDERSTANDING THE UPPER EXTREMITY

