

SLiK Fix is a screw-in tenodesis system recommended for use in both large and small-joint repairs. SLiK Fix is ideally suited for biceps repair.

Features & Benefits

Activation button

- Allows for one-handed insertion which saves time in the operating room

Suture retention tip and tendon washer

- Effective manipulation of tendon

Designed to advance tendon prior to fixation

- Compression fixation along full length of implant



SLiK Fix PEEK CF Screws

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10995	7	15	SLiK Fix, tendonesis screw	PEEK CF	1
10996	8	15	SLiK Fix, tendonesis screw	PEEK CF	1
10997	9	15	SLiK Fix, tendonesis screw	PEEK CF	1

Instrumentation for SLiK Fix Screws

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
11016U	7	Shoulder	SLiK Fix, drill bit	Unsterile	Single-use
11017U	8	Shoulder	SLiK Fix, drill bit	Unsterile	Single-use
11018U	9	Shoulder	SLiK Fix, drill bit	Unsterile	Single-use
11019U	10	Shoulder	SLiK Fix, drill bit	Unsterile	Single-use



Parcus Medical has joined **Anika**



Parcus Medical, LLC.

6423 Parkland Dr, Sarasota, FL

1-941-755-7965 · www.parcusmedical.com

ParcusCS@anika.com