

Twist Knotless

Knotless Screw-in Suture Anchors

Twist Knotless was designed to be extremely reliable and strong. The fully threaded screw-in knotless anchor is ideally suited for rotator cuff repair and Achilles repair.

Features & Benefits

Knotless

- Provides a step-saving alternative to conventional "knotted" suture anchors
- Eliminates "knot stacks" associated with soft tissue irritation

Fully-threaded design

- Optimal cortical fixation
- Reduce the risk of "pull-back"

Make of carbon fiber-reinforced PEEK-OPTIMA

- Non-absorbable, radiolucent, MR safe



Twist Knotless					
Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
11165	4.75	15	Twist Knotless, screw-in, suture anchor, w/suture passer	PEEK	1
11166	5.5	15	Twist Knotless, screw-in, suture anchor, w/suture passer	PEEK	1
11167	6.25	15	Twist Knotless, screw-in, suture anchor, w/suture passer	PEEK	1

Twist Knotless DEX					
Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
11259	4.75	15	Twist Knotless DEX, screw-in, suture anchor, w/suture passer	PEEK	1
11260	5.5	15	Twist Knotless DEX, screw-in, suture anchor, w/suture passer	PEEK	1

Twist Knotless Instrumentation					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
11172	-	Shoulder	Awl, w/taper from 3.6mm to 4.1mm	Unsterile	Reusable
11000U	4	Small Joint	Drill bit	Unsterile	Reusable
11241	4	Small Joint	4.0 x 100mm drill guide, 6-point	Unsterile	Reusable



Parcus Medical has joined **Anika**



Parcus Medical, LLC.

6423 Parkland Dr, Sarasota, FL

1-941-755-7965 · www.parcusmedical.com

ParcusCS@anika.com